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Use a special phosphate bounded investment free of graphite (type FUSION) - Heat the ring at 850°c—900° c—1562 ° F 1652 ° F—Maintain it to this température during 40 to 45 min according to the size of the ring , cast the piece. Let it air cool after casting .

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Melting of ALCO63 is performed in a clean crucible—It can be achieved by high or medium frequency induction or by oxyacetylene or oxypropane flame. In case of melting by flame , the flame has to be adjusted to a neutral and a non recarburating state.

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Choose the size of ingots according to the size of the frame for using only necessary alloy— Switch on the melting machine, as soon as the ingots are wholly melted , start centrifugal action.

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Light the flame and adjust the pressure in order to get a blue point about 4mm long . Heat the ingots in giving to the flame a circular move along the crucible and start the centrifugal action. - Let the ring air cool and put away the piece.

PREPARATION DE LA PIECE - PREPARING OF THE PART

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Proceed as per the conventional way, by scraping and grinding the asperities of the part . Sand blast under high pressure on which the porcelain is to be applied with an abrasive aluminium 100µ—120µ . Ultrasonically clean with distilled water during 5—10 min or employ jet of steam during 5 min to clean the frame surface .

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It is not necessary to make oxydation nor degazing—Opaques are at that time applied on the metallic heat surface according to the recommendations of ceramic manufacturer - ALCO63 has been tested positively on main ceramic on the market.

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Use a special phosphate bounded investment free of graphite (type POMOS ER) - Heat the ring at 900°c—1650 ° c F—Maintain it to this température during 45 min at least according to the size of The ring and cast the piece

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The process is the same as for the realisation of metallo-ceramic prosthesis.

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Break the ring , sand the piece , scrap it with a black diamond bur for beginning , red and then proceed to the finition with a diamond bur yellow ring . Make then le polishing with a rubber wheel and the shining with a felt or a polishing brush and polishing past (Brilliance)

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Ultima garanties optimum propriétés during initial melting of its virgin dental alloys . Ultima Declines all responsibility in the following cases :

- . Remelting of ingots for recycling by dental laboratories
- . Mixing of several heats or grades by the dental laboratories when preparing a melting charge.

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Cobalt chrome alloys seldom cause dermatitis on sensitive subjects, it is recommended to make a contact test .

Verify before the application of the prosthesis if there are other implants in the mouth of the patient . If different alloys are present a "pile effect" may occur .

Our supply are identified by a batch number which must be written on the patient file according tracability process - When casting and polishing , use an adequate vacuum system— technicians are invited to wear a dust mask.