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PROPRIETES MECANIKES *Mechanical properties*

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Use a special phosphate bounded investment free of graphite (type FUSION) - Heat the ring at 850°c—1562 ° F— Maintain it to this température during 40 to 45 min according to the size of the ring , cast the piece. Let it air cool après la coulée .

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Melting of NCU is performed in a clean crucible—It can be achieved by high or medium frequency induction or by oxyacetylene or oxypropane flame. In case of melting by flame , the flame has to be adjusted to a neutral and a non recarburating state.

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Choose the number of ingots according to the size of the frame for using only necessary alloy— Switch on the melting machine, as soon as the ingots are wholly melted , start centrifugal action.

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Light the flame and adjust the pressure in order to get a blue point about 4mm long . Heat the ingots in giving to the flame a circular move along the crucible and start the centrifugal action. - Let the ring air cool and put away the piece.

PREPARATION DE LA PIECE - PREPARING OF THE PART

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It is not necessary to make oxydation nor degazing—Opaques are at that time applied on the metallic surface according to the recommendations of ceramic manufacturer - NCU has been tested positively on main ceramic on the market.

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Ultima guarantees optimum properties during initial melting of its virgin dental alloys . Ultima declines all responsibility in the following cases :

. Remelting of ingots for recycling by dental laboratories

. Mixing of several heats or grades by the dental laboratories when preparing a melting charge.

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This alloy must not be used for patients who are allergic to Nickel .

Verify before the application of the prosthesis if there are other implants in the mouth of the patient . If different alloys are present a "pile effect" may occur .

Our supply are identified by a batch number which must be written on the patient file according to the traceability process .

When casting and polishing use an adequate vacuum system— technicians are invited to wear a dust mask.

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