

Modelling waxes

PROTAX

OSTRALE
DENTAL PRODUCTS



Description

Modelling (base plate) wax sheets are used by the technician to create the preliminary base for the dentist to try in the patient's mouth .

It is recommended for taking occlusion—It may be also used to position the artificial teeth on the metal framework

Features

Color : Peach

Sheet size : 189 mm x 89 mm x 1.5mm

**Pack size : 500 g (24 sheets)
2.5 kg**

Prime softening temperature : 39° - 42° c

Type 1—sot classification

- will not irritate oral tissues
- No unpleasant flavour

Benefits

- Softens without becoming flaky or crumbly and without adhering to fingers
- Trims easily and cleanly
- Shows a smooth glossy surface after gentle flaming
- Does not distort unduly in the mouth
- Easy to soften and manipulate
- No residue on ceramic or synthetic polymer teeth
- Colouring agents will not separate from wax or impregnate plaster

Distribué par Ultima

ZA La Sazerolle—49140 Seiches sur le Loir -
Tél (33)2.41.18.09.88 - Fax (33)2.41.18.09.87
Email : ultima.lp@wanadoo.fr

PROT-500 - 500 g
PROT-2.5—2.5 kg